

Product Datasheet

PKC theta Antibody (453416) [PE] FAB4368P

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB4368P

Updated 8/3/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB4368P



FAB4368P

PKC theta Antibody (453416) [PE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	453416
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	PE
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	5588
Gene Symbol	PRKCQ
Species	Human
Specificity/Sensitivity	Detects human PKCθ. This antibody also detects recombinant human PKCθ, but does not detect any other recombinant human PKC family member.
Immunogen	<i>E. coli</i> -derived recombinant human PKCθ His580-Ser706 Accession # Q04759

Product Application Details	
Applications	Western Blot, Simple Western, CyTOF-ready, Intracellular Staining by Flow Cytometry
Recommended Dilutions	Western Blot, Simple Western, Intracellular Staining by Flow Cytometry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB4368P

H00005588-Q01-10ug	Recombinant Human PKC theta GST (N-Term) Protein
202-IL-010	IL-2 [Unconjugated]
NBP2-10420	PKC theta Overexpression Lysate
NBP1-19371	CD4 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB4368P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

