

# Product Datasheet

## GPR56 Antibody (776123) [PerCP] FAB4636C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB4636C](http://www.novusbio.com/FAB4636C)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB4636C](http://www.novusbio.com/reviews/destination/FAB4636C)



**FAB4636C**

GPR56 Antibody (776123) [PerCP]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.           |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Monoclonal  |
| Clone                       | 776123  |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG2b   |
| Conjugate                   | PerCP   |
| Purity                      | Protein A or G purified from hybridoma culture supernatant  |
| Buffer                      | PBS   |
| Product Description         |   |
| Host                        | Mouse   |
| Gene ID                     | 9289  |
| Species                     | Human   |
| Specificity/Sensitivity     | Detects human GPR56 in direct ELISA.  |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human GPR56 Met1-Val342<br>Accession # Q9Y653 |
| Product Application Details |   |
| Applications                | Flow Cytometry, CyTOF-ready   |
| Recommended Dilutions       | Flow Cytometry, CyTOF-ready   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB4636C**

---

|               |  |
|---------------|--|
| NBP2-27231PCP | Mouse IgG2b Isotype Control (MPC-11) [PerCP] |
| DVE00         | VEGF [HRP]                                   |
| 4634-GP-050   | GPR56  |
| AF3820        | Reelin Antibody [Unconjugated]               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB4636C](http://www.novusbio.com/reviews/submit/FAB4636C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

