

# Product Datasheet

## **RARRES3 Antibody (659209) [Janelia Fluor® 525] FAB6139JF525**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6139JF525](http://www.novusbio.com/FAB6139JF525)

Updated 5/10/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6139JF525](http://www.novusbio.com/reviews/destination/FAB6139JF525)



**FAB6139JF525**

RARRES3 Antibody (659209) [Janelia Fluor® 525]

| Product Information         |  |
|-----------------------------|--|
| Unit Size                   | 0.1 ml   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Monoclonal   |
| Clone                       | 659209   |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG2b  |
| Conjugate                   | Janelia Fluor 525  |
| Purity                      | Protein A or G purified from hybridoma culture supernatant   |
| Buffer                      | 50mM Sodium Borate   |
| Product Description         |  |
| Host                        | Mouse  |
| Gene ID                     | 5920   |
| Gene Symbol                 | RARRES3  |
| Species                     | Human  |
| Specificity/Sensitivity     | Recognizes human RARRES3 in Western blots and direct ELISAs. In Western blots and direct ELISAs, no cross-reactivity with recombinant human (rh) RARRES1 or rhLatexin is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant human RARRES3<br>Asp12-Lys134<br>Accession # Q9UL19  |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.  |
| Product Application Details |  |
| Applications                | Immunocytochemistry  |
| Recommended Dilutions       | Immunocytochemistry  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB6139JF525**

---

|            |  |
|------------|--|
| 210-TA-005 | TNF-alpha [Unconjugated]               |
| NBP3-40346 | Human RARRES3 ELISA Kit (Colorimetric) |
| NBP2-24875 | TLR3 Antibody (40C1285.6) - BSA Free   |
| NBP2-24906 | TLR7 Antibody - BSA Free               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6139JF525](http://www.novusbio.com/reviews/submit/FAB6139JF525)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

