

Product Datasheet

MS4A3 Antibody (489433) [mFluor Violet 610 SE] FAB6247MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/FAB6247MFV610

Updated 8/3/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/FAB6247MFV610



FAB6247MFV610

MS4A3 Antibody (489433) [mFluor Violet 610 SE]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 489433 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | mFluor Violet 610 SE |
| Purity | Protein A or G purified from hybridoma culture supernatant |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|---|
| Host | Mouse |
| Gene ID | 932 |
| Gene Symbol | MS4A3 |
| Species | Human |
| Specificity/Sensitivity | Detects human MS4A3 in direct ELISAs. |
| Immunogen | NS0 mouse myeloma cell line transfected with human MS4A3 Accession # NP_006129 |
| Notes | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Flow Cytometry, CyTOF-ready |
| Recommended Dilutions | Flow Cytometry, CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to FAB6247MFV610

| | |
|--------------------|---|
| H00000932-P02-10ug | Recombinant Human MS4A3 GST (N-Term) Protein |
| H00245802-P01-10ug | Recombinant Human MS4A6E GST (N-Term) Protein |
| H00000932-P01-2ug | Recombinant Human MS4A3 GST (N-Term) Protein |
| AF4117 | CD34 Antibody [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/FAB6247MFV610

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

