

# Product Datasheet

## **EDNRA/Endothelin R Type A Antibody (485711) [DyLight 650] FAB6538W**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6538W](http://www.novusbio.com/FAB6538W)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6538W](http://www.novusbio.com/reviews/destination/FAB6538W)



**FAB6538W**

EDNRA/Endothelin R Type A Antibody (485711) [DyLight 650]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 485711  |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG2b   |
| Conjugate           | DyLight 650   |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate  |

| Product Description     |   |
|-------------------------|---|
| Host                    | Mouse   |
| Gene ID                 | 1909  |
| Gene Symbol             | EDNRA   |
| Species                 | Human   |
| Specificity/Sensitivity | Detects human EDRNA in direct ELISAs.   |
| Immunogen               | NS0 mouse myeloma cell line transfected with human EDRNA<br>Accession # NP_001948 |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Immunohistochemistry   |
| Recommended Dilutions       | Immunohistochemistry   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB6538W**

|                   |  |
|-------------------|--|
| NBP2-27231C       | Mouse IgG2b Isotype Control (MPC-11) [DyLight 650]               |
| H00001909-P01-2ug | Recombinant Human EDNRA/Endothelin R Type A GST (N-Term) Protein |
| 210-TA-005        | TNF-alpha [Unconjugated]   |
| NBP2-80131        | Human EDNRA/Endothelin R Type A ELISA Kit (Colorimetric)         |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6538W](http://www.novusbio.com/reviews/submit/FAB6538W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

