

# Product Datasheet

## **CCL3/MIP-1 alpha Antibody (755012) [Allophycocyanin] FAB6625A**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB6625A](http://www.novusbio.com/FAB6625A)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB6625A](http://www.novusbio.com/reviews/destination/FAB6625A)



**FAB6625A**

CCL3/MIP-1 alpha Antibody (755012) [Allophycocyanin]

| Product Information         |                                                                                                                                                                                                                                                                                                                        |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                               |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                             |
| Clone                       | 755012                                                                                                                                                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                     |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                                                                                  |
| Conjugate                   | Allophycocyanin                                                                                                                                                                                                                                                                                                        |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                                                                                                                                                                                                                             |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                    |
| Product Description         |                                                                                                                                                                                                                                                                                                                        |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                  |
| Gene ID                     | 6348                                                                                                                                                                                                                                                                                                                   |
| Gene Symbol                 | CCL3                                                                                                                                                                                                                                                                                                                   |
| Species                     | Rat                                                                                                                                                                                                                                                                                                                    |
| Specificity/Sensitivity     | Detects rat CCL3/MIP-1 $\alpha$ in direct ELISAs and Western blots. In direct ELISAs, less than 20% cross-reactivity with recombinant mouse CCL3 and recombinant cotton rat CCL3 is observed and no cross-reactivity with recombinant human CCL3, recombinant canine CCL3, or recombinant rat CCL4/MIP-1b is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant rat CCL3/MIP-1 $\alpha$<br>Ala24-Ala92<br>Accession # P50229                                                                                                                                                                                                                       |
| Product Application Details |                                                                                                                                                                                                                                                                                                                        |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                           |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                 |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB6625A**

---

|                  |                                                        |
|------------------|--------------------------------------------------------|
| IC0041A          | Mouse IgG2b Isotype Control (133303) [Allophycocyanin] |
| NBP2-35193-100ug | Recombinant Mouse CCL3/MIP-1 alpha Protein             |
| 210-TA-005       | TNF-alpha [Unconjugated]                               |
| MMA00            | CCL3/MIP-1 alpha [HRP]                                 |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB6625A](http://www.novusbio.com/reviews/submit/FAB6625A)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

