

# Product Datasheet

## MS4A4A Antibody (818112) [Janelia Fluor® 646] FAB7797J

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB7797J](http://www.novusbio.com/FAB7797J)

Updated 7/29/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB7797J](http://www.novusbio.com/reviews/destination/FAB7797J)



**FAB7797J**

MS4A4A Antibody (818112) [Janelia Fluor® 646]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 818112                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2b                                                                                       |
| Conjugate           | Janelia Fluor 646                                                                           |
| Purity              | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                       |
|-------------------------|---------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                 |
| Gene ID                 | 51338                                                                                 |
| Gene Symbol             | MS4A4A                                                                                |
| Species                 | Human                                                                                 |
| Specificity/Sensitivity | Detects human MS4A4A.                                                                 |
| Immunogen               | Mouse myeloma cell line NS0 transfected with human MS4A4A<br>Accession # Q96JQ5       |
| Notes                   | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus. |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Flow Cytometry, CyTOF-ready                                            |
| Recommended Dilutions       | Flow Cytometry, CyTOF-ready                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to FAB7797J**

---

|                   |                                               |
|-------------------|-----------------------------------------------|
| H00051338-P01-2ug | Recombinant Human MS4A4A GST (N-Term) Protein |
| 7268-CT-100       | CTLA-4 [Unconjugated]                         |
| NBL1-13318        | MS4A4A Overexpression Lysate                  |
| NBP2-22203        | ERK1 Antibody (1E5) - BSA Free                |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB7797J](http://www.novusbio.com/reviews/submit/FAB7797J)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

