

# Product Datasheet

## **HHEX Antibody (2018B) [DyLight 650] FAB83771W**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB83771W](http://www.novusbio.com/FAB83771W)

Updated 10/7/2024 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB83771W](http://www.novusbio.com/reviews/destination/FAB83771W)



**FAB83771W**

HHEX Antibody (2018B) [DyLight 650]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | 2018B                                                                                       |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG                                                                                         |
| Conjugate                   | DyLight 650                                                                                 |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                  |
| Buffer                      | 50mM Sodium Borate                                                                          |
| Product Description         |                                                                                             |
| Host                        | Rabbit                                                                                      |
| Gene ID                     | 3087                                                                                        |
| Species                     | Human, Mouse, Rat                                                                           |
| Specificity/Sensitivity     | Detects human HHEX in direct ELISAs and Western blots.                                      |
| Immunogen                   | <i>E. coli</i> -derived recombinant human HHEX<br>Thr111-Gly270<br>Accession # Q03014       |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.           |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry                                           |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to FAB83771W**

---

|             |                                          |
|-------------|------------------------------------------|
| NBP2-24891C | Rabbit IgG Isotype Control [DyLight 650] |
| NBP2-23108  | Recombinant Human HHEX His Protein       |
| 2496-ZN-010 | Insulysin/IDE [Unconjugated]             |
| NBL1-11531  | HHEX Overexpression Lysate               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB83771W](http://www.novusbio.com/reviews/submit/FAB83771W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

