

# Product Datasheet

## Glypican 4 Antibody (961609) [HRP] FAB9195H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB9195H](http://www.novusbio.com/FAB9195H)

Updated 10/7/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB9195H](http://www.novusbio.com/reviews/destination/FAB9195H)



**FAB9195H**

Glypican 4 Antibody (961609) [HRP]

| Product Information         |                                                                                                            |
|-----------------------------|------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                |
| Storage                     | Store at 4C in the dark.                                                                                   |
| Clonality                   | Monoclonal                                                                                                 |
| Clone                       | 961609                                                                                                     |
| Preservative                | No Preservative                                                                                            |
| Isotype                     | IgG1                                                                                                       |
| Conjugate                   | HRP                                                                                                        |
| Purity                      | Protein A or G purified from hybridoma culture supernatant                                                 |
| Buffer                      | PBS                                                                                                        |
| Product Description         |                                                                                                            |
| Host                        | Mouse                                                                                                      |
| Gene ID                     | 2239                                                                                                       |
| Species                     | Human                                                                                                      |
| Specificity/Sensitivity     | Detects human Glypican 4 in direct ELISAs.                                                                 |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human Glypican 4 Met1-Ser529<br>Accession # O75487 |
| Product Application Details |                                                                                                            |
| Applications                | ELISA, Flow Cytometry, CyTOF-ready                                                                         |
| Recommended Dilutions       | Flow Cytometry, ELISA, CyTOF-ready                                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB9195H**

---

|               |                                        |
|---------------|----------------------------------------|
| NBP1-97074    | Mouse IgG1 Isotype Control (MG1) [HRP] |
| NBP3-17858PEP | Glypican 4 Recombinant Protein Antigen |
| 233-FB-025    | FGF basic/FGF2/bFGF [Unconjugated]     |
| 9195-GP-050   | Glypican 4 [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB9195H](http://www.novusbio.com/reviews/submit/FAB9195H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

