

# Product Datasheet

## ABCG2/CD338 Antibody (5D3) [HRP] FAB995H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/FAB995H](http://www.novusbio.com/FAB995H)

Updated 8/1/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/FAB995H](http://www.novusbio.com/reviews/destination/FAB995H)



**FAB995H**

ABCG2/CD338 Antibody (5D3) [HRP]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Monoclonal  |
| Clone                       | 5D3   |
| Preservative                | No Preservative   |
| Isotype                     | IgG2b   |
| Conjugate                   | HRP   |
| Purity                      | 0   |
| Buffer                      | PBS   |
| Product Description         |   |
| Host                        | Mouse   |
| Gene ID                     | 9429  |
| Gene Symbol                 | ABCG2   |
| Species                     | Human   |
| Specificity/Sensitivity     | Detects human ABCG2 in flow cytometry and immunocytochemistry.                              |
| Immunogen                   | 3T3 cells transduced with human ABCG2   |
| Product Application Details |   |
| Applications                | Flow Cytometry, CyTOF-ready, Immunocytochemistry  |
| Recommended Dilutions       | Flow Cytometry, Immunocytochemistry, CyTOF-ready  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to FAB995H**

---

|                   |  |
|-------------------|--|
| NBP2-27231H       | Mouse IgG2b Isotype Control (MPC-11) [HRP]         |
| H00009429-P01-2ug | Recombinant Human ABCG2/CD338 GST (N-Term) Protein |
| NBP3-27492        | Human ABCG2/CD338 ELISA Kit (Colorimetric)         |
| AF231             | EGFR Antibody [Unconjugated]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/FAB995H](http://www.novusbio.com/reviews/submit/FAB995H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

