

# Product Datasheet

## SAH3 Antibody [Janelia Fluor® 635] NB100-400JF635

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-400JF635](http://www.novusbio.com/NB100-400JF635)

Updated 8/20/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-400JF635](http://www.novusbio.com/reviews/destination/NB100-400JF635)



**NB100-400JF635**

SAH3 Antibody [Janelia Fluor® 635]

| Product Information         |                                                                                                                                                                                                  |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                      |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                         |
| Clonality                   | Polyclonal                                                                                                                                                                                       |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                                              |
| Conjugate                   | Janelia Fluor 635                                                                                                                                                                                |
| Purity                      | Immunogen affinity purified                                                                                                                                                                      |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                               |
| Product Description         |                                                                                                                                                                                                  |
| Host                        | Rabbit                                                                                                                                                                                           |
| Gene ID                     | 23382                                                                                                                                                                                            |
| Gene Symbol                 | AHCYL2                                                                                                                                                                                           |
| Species                     | Human                                                                                                                                                                                            |
| Immunogen                   | The immunogen recognized by this antibody maps to a region between residues 1 and 50 of human adenosylhomocysteinase 3 (KIAA0828) using the numbering given in entry NP_056143.1 (GeneID 23382). |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                                            |
| Product Application Details |                                                                                                                                                                                                  |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                                |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                                                                                                                                |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                           |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NB100-400JF635**

---

|                    |                                             |
|--------------------|---------------------------------------------|
| H00023382-P01-10ug | Recombinant Human SAH3 GST (N-Term) Protein |
| 6507-IL-010/CF     | IL-4 [Unconjugated]                         |
| NBP2-47293PEP      | SAH3 Recombinant Protein Antigen            |
| 7954-GM-010/CF     | GM-CSF [Unconjugated]                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-400JF635](http://www.novusbio.com/reviews/submit/NB100-400JF635)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

