

# Product Datasheet

## Hepatitis B preS2 Antibody [DyLight 405] NB100-93580V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-93580V](http://www.novusbio.com/NB100-93580V)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-93580V](http://www.novusbio.com/reviews/destination/NB100-93580V)



**NB100-93580V**

Hepatitis B preS2 Antibody [DyLight 405]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                           |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Polyclonal  |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG   |
| Conjugate                   | DyLight 405   |
| Purity                      | Immunogen affinity purified   |
| Buffer                      | 50mM Sodium Borate  |
| Product Description         |   |
| Host                        | Rabbit  |
| Species                     | Human   |
| Reactivity Notes            | 0   |
| Specificity/Sensitivity     | Hepatitis B Virus preS2 antigen   |
| Immunogen                   | Synthetic peptide within residues 130-150 of the 42 kDa preS2 HBV protein. This sequence represents a B cell epitope. |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                     |
| Product Application Details |   |
| Applications                | Western Blot  |
| Recommended Dilutions       | Western Blot  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NB100-93580V**

---

|                 |  |
|-----------------|--|
| NBP2-24891V     | Rabbit IgG Isotype Control [DyLight 405] |
| NBP2-35884-10ug | Recombinant Hepatitis B preS2 Protein    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-93580V](http://www.novusbio.com/reviews/submit/NB100-93580V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

