

# Product Datasheet

## Goat anti-Guinea Pig IgG (H+L) Secondary Antibody [Allophycocyanin] NB120-7134APC

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB120-7134APC](http://www.novusbio.com/NB120-7134APC)

Updated 9/6/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB120-7134APC](http://www.novusbio.com/reviews/destination/NB120-7134APC)



**NB120-7134APC**

Goat anti-Guinea Pig IgG (H+L) Secondary Antibody [Allophycocyanin]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                      |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                                                                               |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                          |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                             |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                           |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                   |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                  |
| Conjugate                   | Allophycocyanin                                                                                                                                                                                                                                                                                                                                      |
| Purity                      | Multi-step                                                                                                                                                                                                                                                                                                                                           |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                                  |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                      |
| Host                        | Goat                                                                                                                                                                                                                                                                                                                                                 |
| Species                     | Guinea Pig                                                                                                                                                                                                                                                                                                                                           |
| Specificity/Sensitivity     | <p>This antibody was pre-adsorbed against Bovine, Chicken, Goat, Hamster, Horse, Human, Mouse, Rabbit, Rat, and Sheep Serum Proteins.</p> <p>No reaction was observed against Bovine, Chicken, Goat, Hamster, Horse, Human, Mouse, Rabbit, Rat and Sheep Serum Proteins. Minimal cross reactivity was observed against the species listed above.</p> |
| Immunogen                   | Guinea Pig IgG whole molecule                                                                                                                                                                                                                                                                                                                        |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                      |
| Applications                | Western Blot, Dot Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                                                                                                                                                                                                                                         |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Dot Blot                                                                                                                                                                                                                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB120-7134APC](http://www.novusbio.com/reviews/submit/NB120-7134APC)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

