

# Product Datasheet

## **Troponin I Type 1 (slow skeletal) Antibody (7G2) [Janelia Fluor® 525] NB200-432JF525**

Unit Size: 0.2 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB200-432JF525](http://www.novusbio.com/NB200-432JF525)

Updated 8/20/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB200-432JF525](http://www.novusbio.com/reviews/destination/NB200-432JF525)



**NB200-432JF525**

Troponin I Type 1 (slow skeletal) Antibody (7G2) [Janelia Fluor® 525]

Product Information	
Unit Size	0.2 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	7G2
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	Janelia Fluor 525
Purity	Protein A purified
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	7135
Gene Symbol	TNNI1
Species	Human, Invertebrate
Reactivity Notes	Invertebrate reactivity reported in scientific literature (PMID: 26213227).
Specificity/Sensitivity	Troponin I from human skeletal muscle. Cross-reaction of this with cardiac troponin I is less than 1.5%.
Immunogen	Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with Human skTnl.
Notes	Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.

Product Application Details	
Applications	Western Blot, Immunoassay, Sandwich ELISA, Immunoaffinity Purification
Recommended Dilutions	Western Blot, Immunoassay, Sandwich ELISA, Immunoaffinity Purification
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB200-432JF525**

---

NBP2-22892	Recombinant Human Troponin I Type 1 (slow skeletal) His Protein
NBP2-76598	Human Troponin I Type 1 (slow skeletal) ELISA Kit (Colorimetric)
MAB4470	Myosin Heavy Chain Antibody (MF20) [Unconjugated]
NB110-68136	Tenascin C Antibody (4C8MS) - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB200-432JF525](http://www.novusbio.com/reviews/submit/NB200-432JF525)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

