

# Product Datasheet

## QKI/Quaking Antibody [HRP] NB300-240H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB300-240H](http://www.novusbio.com/NB300-240H)

Updated 12/12/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB300-240H](http://www.novusbio.com/reviews/destination/NB300-240H)



**NB300-240H**

QKI/Quaking Antibody [HRP]

| Product Information         |                                                                                                                                                                                                                                                                                                                     |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                         |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                            |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                          |
| Preservative                | No Preservative                                                                                                                                                                                                                                                                                                     |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                 |
| Conjugate                   | HRP                                                                                                                                                                                                                                                                                                                 |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                         |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                 |
| Product Description         |                                                                                                                                                                                                                                                                                                                     |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                              |
| Gene ID                     | 9444                                                                                                                                                                                                                                                                                                                |
| Gene Symbol                 | QKI                                                                                                                                                                                                                                                                                                                 |
| Species                     | Human, Mouse                                                                                                                                                                                                                                                                                                        |
| Immunogen                   | The immunogen this antibody was made to, maps to the C-terminus of isoform 5 of human quaking homolog according to entry AAH19917.1 (GeneID 9444). The epitope is not found in isoforms 1, 2, 3, 4, 6, 7 or 7B according to TrEMBL entries Q9P0Y1, Q9P0Y0, Q9P0X6, Q9P0X8, Q96PU7, Q969L9 and Q96PU6, respectively. |
| Product Application Details |                                                                                                                                                                                                                                                                                                                     |
| Applications                | Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation                                                                                                                                                                                                                                            |
| Recommended Dilutions       | Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                              |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB300-240H**

---

|             |                                           |
|-------------|-------------------------------------------|
| NBP2-24891H | Rabbit IgG Isotype Control [HRP]          |
| NBP2-23390  | Recombinant Human QKI/Quaking His Protein |
| 538-MG-100  | MAG/Siglec-4a                             |
| NBL1-15028  | QKI/Quaking Overexpression Lysate         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB300-240H](http://www.novusbio.com/reviews/submit/NB300-240H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

