

# Product Datasheet

## Goat anti-Human IgG (H+L) Secondary Antibody [DyLight 350] NB7483UV

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB7483UV](http://www.novusbio.com/NB7483UV)

Updated 9/6/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB7483UV](http://www.novusbio.com/reviews/destination/NB7483UV)



**NB7483UV**

Goat anti-Human IgG (H+L) Secondary Antibody [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Purity                      | Immunogen affinity purified                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Host                        | Goat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Specificity/Sensitivity     | By immunoelectrophoresis and ELISA this Goat anti-Human IgG (H+L) Secondary Antibody (Pre-adsorbed) reacts specifically with human IgG and with light chains common to other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Antiserum was cross adsorbed using bovine, chicken, horse, mouse, pig, rabbit and rat immunosorbents to remove cross reactive Antibodies. Less than 0.1% cross reactivity to bovine, chicken, horse, mouse, pig, rabbit and rat IgG was detected. This antibody may cross react with IgG from other species. |
| Immunogen                   | This Goat anti-Human IgG (H+L) Secondary Antibody (Pre-adsorbed) was developed against human IgG.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Applications                | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB7483UV](http://www.novusbio.com/reviews/submit/NB7483UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

