

# Product Datasheet

## FKBP25 Antibody [FITC] NBP1-41071F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-41071F](http://www.novusbio.com/NBP1-41071F)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-41071F](http://www.novusbio.com/reviews/destination/NBP1-41071F)



**NBP1-41071F**

FKBP25 Antibody [FITC]

| Product Information         |                                                                                                                                                                           |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                               |
| Storage                     | Store at 4C in the dark.                                                                                                                                                  |
| Clonality                   | Polyclonal                                                                                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                        |
| Isotype                     | IgG                                                                                                                                                                       |
| Conjugate                   | FITC                                                                                                                                                                      |
| Purity                      | Immunogen affinity purified                                                                                                                                               |
| Buffer                      | PBS                                                                                                                                                                       |
| Product Description         |                                                                                                                                                                           |
| Host                        | Rabbit                                                                                                                                                                    |
| Gene ID                     | 2287                                                                                                                                                                      |
| Gene Symbol                 | FKBP3                                                                                                                                                                     |
| Species                     | Human                                                                                                                                                                     |
| Immunogen                   | The immunogen for this product maps to a region between residue 50 and 100 of human FK506 binding protein 3 using the numbering given in entry NP_002004.1 (GeneID 2287). |
| Product Application Details |                                                                                                                                                                           |
| Applications                | Immunoprecipitation, Western Blot (Negative)                                                                                                                              |
| Recommended Dilutions       | Immunoprecipitation, Western Blot (Negative)                                                                                                                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP1-41071F**

---

|            |                                         |
|------------|-----------------------------------------|
| NBP2-24892 | Rabbit IgG Isotype Control [FITC]       |
| NBP1-41160 | Recombinant Human FKBP25 Protein        |
| NBL1-10737 | FKBP25 Overexpression Lysate            |
| NB600-1397 | Transportin 1 Antibody (D45) - BSA Free |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-41071F](http://www.novusbio.com/reviews/submit/NBP1-41071F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

