

# Product Datasheet

## LAMP-1/CD107a Antibody [Allophycocyanin] NBP1-77241APC

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-77241APC](http://www.novusbio.com/NBP1-77241APC)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-77241APC](http://www.novusbio.com/reviews/destination/NBP1-77241APC)



**NBP1-77241APC**

LAMP-1/CD107a Antibody [Allophycocyanin]

| Product Information         |                                                                                                                                                                                                           |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                               |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                  |
| Clonality                   | Polyclonal                                                                                                                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                        |
| Isotype                     | IgG                                                                                                                                                                                                       |
| Conjugate                   | Allophycocyanin                                                                                                                                                                                           |
| Purity                      | Peptide affinity purified                                                                                                                                                                                 |
| Buffer                      | PBS                                                                                                                                                                                                       |
| Product Description         |                                                                                                                                                                                                           |
| Host                        | Rabbit                                                                                                                                                                                                    |
| Gene ID                     | 3916                                                                                                                                                                                                      |
| Gene Symbol                 | LAMP1                                                                                                                                                                                                     |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                         |
| Marker                      | Late Endosome Marker                                                                                                                                                                                      |
| Specificity/Sensitivity     | LAMP-1 antibody is predicted to not cross-react with LAMP-2. Due to extensive post-translational modification, LAMP-1 often migrates at higher than predicted molecular weight in SDS-PAGE.               |
| Immunogen                   | Antibody was raised against a 15 amino acid synthetic peptide from near the center of human LAMP-1. The immunogen is located within amino acids 230 - 280 of LAMP-1. Amino Acid Sequence: NLTYERKDNTTVTRL |
| Product Application Details |                                                                                                                                                                                                           |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen                                                             |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP1-77241APC**

---

|               |                                                       |
|---------------|-------------------------------------------------------|
| IC1051A       | Rabbit IgG Isotype Control (60024B) [Allophycocyanin] |
| NBP2-61619PEP | LAMP-1/CD107a Recombinant Protein Antigen             |
| 210-TA-005    | TNF-alpha [Unconjugated]                              |
| 4320-LM-050   | LAMP-1/CD107a                                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-77241APC](http://www.novusbio.com/reviews/submit/NBP1-77241APC)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

