

# Product Datasheet

## BAIAP3 Antibody [PE] NBP2-41080PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-41080PE](http://www.novusbio.com/NBP2-41080PE)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-41080PE](http://www.novusbio.com/reviews/destination/NBP2-41080PE)



**NBP2-41080PE**

BAIAP3 Antibody [PE]

| Product Information         |                                                                                                                                                                                                           |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                               |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                  |
| Clonality                   | Polyclonal                                                                                                                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                        |
| Isotype                     | IgG                                                                                                                                                                                                       |
| Conjugate                   | PE                                                                                                                                                                                                        |
| Purity                      | Peptide affinity purified                                                                                                                                                                                 |
| Buffer                      | PBS                                                                                                                                                                                                       |
| Product Description         |                                                                                                                                                                                                           |
| Host                        | Rabbit                                                                                                                                                                                                    |
| Gene ID                     | 8938                                                                                                                                                                                                      |
| Gene Symbol                 | BAIAP3                                                                                                                                                                                                    |
| Species                     | Human                                                                                                                                                                                                     |
| Specificity/Sensitivity     | BAP3 is known to have two isoforms; this BAP3 antibody will recognize only isoform 2. Lower molecular weight bands may represent cleavage or degradation products.                                        |
| Immunogen                   | Antibody was raised against a 13 amino acid synthetic peptide from near the amino terminus of human BAP3. The immunogen is located within amino acids 60 - 110 of BAP3. Amino Acid Sequence: CRQSPHPGHTQV |
| Product Application Details |                                                                                                                                                                                                           |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry                                                                                                                                  |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                    |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-41080PE**

---

|               |                                    |
|---------------|------------------------------------|
| NBP2-24983    | Rabbit IgG Isotype Control [PE]    |
| NBP1-88709PEP | BAIAP3 Recombinant Protein Antigen |
| 923-AN-025    | Angiopoietin-1 [Unconjugated]      |
| NBP2-07934    | BAIAP3 Overexpression Lysate       |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-41080PE](http://www.novusbio.com/reviews/submit/NBP2-41080PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

