

# Product Datasheet

## Chromogranin A Antibody (CGA/413 + CHGA/777 + CHGA/798) [HRP] NBP2-47846H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-47846H](http://www.novusbio.com/NBP2-47846H)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-47846H](http://www.novusbio.com/reviews/destination/NBP2-47846H)



**NBP2-47846H**

Chromogranin A Antibody (CGA/413 + CHGA/777 + CHGA/798) [HRP]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | CGA/413 + CHGA/777 + CHGA/798                                                               |
| Preservative                | No Preservative                                                                             |
| Isotype                     | IgG2b Kappa                                                                                 |
| Conjugate                   | HRP                                                                                         |
| Purity                      | Protein A or G purified                                                                     |
| Buffer                      | PBS                                                                                         |
| Product Description         |                                                                                             |
| Host                        | Mouse                                                                                       |
| Gene ID                     | 1113                                                                                        |
| Gene Symbol                 | CHGA                                                                                        |
| Species                     | Human, Mouse, Rat, Porcine, Monkey                                                          |
| Marker                      | Neuroendocrine Marker                                                                       |
| Immunogen                   | Recombinant human Chromogranin A protein (Uniprot: P10645)                                  |
| Product Application Details |                                                                                             |
| Applications                | Immunohistochemistry, Immunohistochemistry-Paraffin                                         |
| Recommended Dilutions       | Immunohistochemistry, Immunohistochemistry-Paraffin                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP2-47846H**

---

|                 |                                                            |
|-----------------|------------------------------------------------------------|
| NBP1-43317H     | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [HRP] |
| NBP2-61453-10ug | Recombinant Human Chromogranin A His Protein               |
| DKK300          | Kallikrein 3/PSA [HRP]                                     |
| DY9098-05       | Chromogranin A [Biotin]                                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-47846H](http://www.novusbio.com/reviews/submit/NBP2-47846H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

