

# Product Datasheet

## Chromogranin A Antibody (CHGA/1731R) [DyLight 594] NBP2-54316DL594

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-54316DL594](http://www.novusbio.com/NBP2-54316DL594)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-54316DL594](http://www.novusbio.com/reviews/destination/NBP2-54316DL594)



**NBP2-54316DL594**

Chromogranin A Antibody (CHGA/1731R) [DyLight 594]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | CHGA/1731R                                                                                  |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Conjugate           | DyLight 594                                                                                 |
| Purity              | Protein A or G purified                                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                   |
|---------------------|-----------------------------------------------------------------------------------|
| Host                | Rabbit                                                                            |
| Gene ID             | 1113                                                                              |
| Gene Symbol         | CHGA                                                                              |
| Species             | Human, Rat (Negative)                                                             |
| Reactivity Notes    | Does not react with Rat.                                                          |
| Marker              | Neuroendocrine Marker                                                             |
| Immunogen           | Recombinant human full-length Chromogranin A protein (Uniprot: P10645)            |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Immunohistochemistry, Immunohistochemistry-Paraffin, Protein Array     |
| Recommended Dilutions       | Immunohistochemistry, Immunohistochemistry-Paraffin, Protein Array     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP2-54316DL594**

---

|                 |                                              |
|-----------------|----------------------------------------------|
| NBP2-24891DL594 | Rabbit IgG Isotype Control [DyLight 594]     |
| NBP2-61453-10ug | Recombinant Human Chromogranin A His Protein |
| DKK300          | Kallikrein 3/PSA [HRP]                       |
| DY9098-05       | Chromogranin A [Biotin]                      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-54316DL594](http://www.novusbio.com/reviews/submit/NBP2-54316DL594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

