

# Product Datasheet

## LIM Kinase 1/2 [p Thr508] Antibody [Alexa Fluor® 532] NBP2-54645AF532

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-54645AF532](http://www.novusbio.com/NBP2-54645AF532)

Updated 12/12/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-54645AF532](http://www.novusbio.com/reviews/destination/NBP2-54645AF532)



**NBP2-54645AF532**

LIM Kinase 1/2 [p Thr508] Antibody [Alexa Fluor® 532]

| Product Information         |                                                                                                           |
|-----------------------------|-----------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.               |
| Storage                     | Store at 4C in the dark.                                                                                  |
| Clonality                   | Polyclonal                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                        |
| Isotype                     | IgG                                                                                                       |
| Conjugate                   | Alexa Fluor 532                                                                                           |
| Purity                      | Epitope affinity purified                                                                                 |
| Buffer                      | 50mM Sodium Borate                                                                                        |
| Product Description         |                                                                                                           |
| Host                        | Rabbit                                                                                                    |
| Gene ID                     | 3984                                                                                                      |
| Gene Symbol                 | LIMK1                                                                                                     |
| Species                     | Human, Mouse, Rat                                                                                         |
| Specificity/Sensitivity     | This antibody detects endogenous level of LIM Kinase 1/2 only when phosphorylated at threonine 508.       |
| Immunogen                   | A peptide sequence around phosphorylation site of threonine 508 (R-Y-T (p) -V-V) derived from Human LIMK1 |
| Product Application Details |                                                                                                           |
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-54645AF532**

---

|                 |                                               |
|-----------------|-----------------------------------------------|
| NBP2-24891AF532 | Rabbit IgG Isotype Control [Alexa Fluor® 532] |
| 236-EG-200      | EGF [Unconjugated]                            |
| AF887           | AKT [p Ser473] Antibody - Pan Specific        |
| NBP2-29463      | p21/CIP1/CDKN1A Antibody (WA-1)               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-54645AF532](http://www.novusbio.com/reviews/submit/NBP2-54645AF532)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

