

# Product Datasheet

## TCF15 Recombinant Protein Antigen NBP2-55818PEP

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-55818PEP](http://www.novusbio.com/NBP2-55818PEP)

Updated 2/4/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-55818PEP](http://www.novusbio.com/reviews/destination/NBP2-55818PEP)



**NBP2-55818PEP****TCF15 Recombinant Protein Antigen**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 100 ul                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Storage</b>                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Preservative</b>                | No Preservative                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>                      | >80% by SDS-PAGE and Coomassie blue staining                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Buffer</b>                      | PBS and 1M Urea, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Description</b>                 | <p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human TCF15</p> <p><b>Source:</b> <i>E.coli</i></p> <p><b>Amino Acid Sequence:</b> CLSNQRKGGGRRDLGGSCCLKVRGVAPLRGPRR</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p>                                                                                                                                     |
| <b>Gene ID</b>                     | 6939                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Gene Symbol</b>                 | TCF15                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Species</b>                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Applications</b>                | Antibody Competition                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Recommended Dilutions</b>       | Antibody Competition 10-100 molar excess                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Application Notes</b>           | <p>This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP2-55818. It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml. For current lot information, including availability, please contact our technical support team click <a href="mailto:nb-technical@bio-techne.com">nb-technical@bio-techne.com</a>.</p> |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-55818PEP**

---

|             |                                 |
|-------------|---------------------------------|
| 423-F8-025  | FGF-8 [Unconjugated]            |
| NBP2-55818  | TCF15 Antibody - BSA Free       |
| NB100-65805 | CD81 Antibody (1D6) - BSA Free  |
| NB100-56511 | MyoD Antibody (5.8A) - BSA Free |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-55818PEP](http://www.novusbio.com/reviews/submit/NBP2-55818PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

