

Product Datasheet

SAP30 Antibody - BSA Free

NBP2-59195-0.05mg

Unit Size: 0.05 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59195

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59195



NBP2-59195-0.05mg

SAP30 Antibody - BSA Free

Product Information

Unit Size	0.05 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide and 0.05% ProClin 300
Isotype	IgG
Purity	Protein G purified
Buffer	PBS

Product Description

Description	Novus Biologicals Rabbit SAP30 Antibody - BSA Free (NBP2-59195) is a polyclonal antibody validated for use in WB and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8819
Gene Symbol	SAP30
Species	Human
Immunogen	The exact immunogen is proprietary information.


Product Application Details

Applications	Western Blot, Chromatin Immunoprecipitation (ChIP), Chromatin Immunoprecipitation Sequencing
Recommended Dilutions	Western Blot 1:1000, Chromatin Immunoprecipitation (ChIP) 2-5 ug, Chromatin Immunoprecipitation Sequencing




Images

Western Blot: SAP30 Antibody [NBP2-59195] - Western blot was performed on nuclear extracts from HeLa cells (20 ug) using the antibody against SAP30 diluted 1:1000 in TBS-Tween containing 5% skimmed milk. The molecular weight marker (in kDa) is shown on the left; the location of the protein of interest is indicated on the right.



Chromatin Immunoprecipitation: SAP30 Antibody [NBP2-59195] - ChIP assays were performed using HeLa cells, the antibody against SAP30 and optimized primer sets for qPCR. ChIP was performed using sheared chromatin from 4 million cells. A titration of the antibody consisting of 1, 2 and 5 ug per ChIP experiment was analysed. IgG (1 ug/ IP) was used as negative IP control. QPCR was performed with primers for the EIF2S3 and BRCA1 promoters, used as positive controls, and for the MYOD1 gene and the Sat2 satellite repeat, used as negative controls. Figure shows the recovery, expressed as a % of input (the relative amount of immunoprecipitated DNA compared to input DNA after qPCR analysis).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59195-0.05mg

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00008819-P01-10ug	Recombinant Human SAP30 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59195

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

