

# Product Datasheet

## OGFOD1 Antibody (OTI1D9) [DyLight 594] NBP2-73145DL594

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-73145DL594](http://www.novusbio.com/NBP2-73145DL594)

Updated 11/9/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-73145DL594](http://www.novusbio.com/reviews/destination/NBP2-73145DL594)



**NBP2-73145DL594**

OGFOD1 Antibody (OTI1D9) [DyLight 594]

| Product Information         |                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.         |
| Storage                     | Store at 4C in the dark.                                                                            |
| Clonality                   | Monoclonal                                                                                          |
| Clone                       | OTI1D9                                                                                              |
| Preservative                | 0.05% Sodium Azide                                                                                  |
| Isotype                     | IgG1                                                                                                |
| Conjugate                   | DyLight 594                                                                                         |
| Purity                      | Immunogen affinity purified                                                                         |
| Buffer                      | 50mM Sodium Borate                                                                                  |
| Product Description         |                                                                                                     |
| Host                        | Mouse                                                                                               |
| Gene ID                     | 55239                                                                                               |
| Gene Symbol                 | OGFOD1                                                                                              |
| Species                     | Human, Monkey                                                                                       |
| Immunogen                   | Full length recombinant protein of human OGFOD1 (NP_060703) produced in HEK293T cells.              |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                   |
| Product Application Details |                                                                                                     |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                              |





### Novus Biologicals USA

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### Bio-Techne Canada

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### Bio-Techne Ltd

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### General Contact Information

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### Products Related to NBP2-73145DL594

|                    |                                                |
|--------------------|------------------------------------------------|
| NBP1-97005DL594    | Mouse IgG1 Isotype Control (MG1) [DyLight 594] |
| H00055239-P01-10ug | Recombinant Human OGFOD1 GST (N-Term) Protein  |
| 210-TA-005         | TNF-alpha [Unconjugated]                       |
| NBL1-13914         | OGFOD1 Overexpression Lysate                   |

### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-73145DL594](http://www.novusbio.com/reviews/submit/NBP2-73145DL594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

