

# Product Datasheet

## t-Plasminogen Activator/tPA Antibody (OTI3H3) [PerCP] NBP2-74568PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-74568PCP](http://www.novusbio.com/NBP2-74568PCP)

Updated 11/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-74568PCP](http://www.novusbio.com/reviews/destination/NBP2-74568PCP)



**NBP2-74568PCP**

t-Plasminogen Activator/tPA Antibody (OTI3H3) [PerCP]

| Product Information         |                                                                                                                      |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                               |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                          |
| Storage                     | Store at 4C in the dark.                                                                                             |
| Clonality                   | Monoclonal                                                                                                           |
| Clone                       | OTI3H3                                                                                                               |
| Preservative                | 0.05% Sodium Azide                                                                                                   |
| Isotype                     | IgG2b                                                                                                                |
| Conjugate                   | PerCP                                                                                                                |
| Purity                      | Immunogen affinity purified                                                                                          |
| Buffer                      | PBS                                                                                                                  |
| Product Description         |                                                                                                                      |
| Host                        | Mouse                                                                                                                |
| Gene ID                     | 5327                                                                                                                 |
| Gene Symbol                 | PLAT                                                                                                                 |
| Species                     | Human                                                                                                                |
| Immunogen                   | Human recombinant protein fragment corresponding to amino acids 157-426 of human PLAT(NP_000921) produced in E.coli. |
| Product Application Details |                                                                                                                      |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                                                    |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-74568PCP**

---

|                 |                                                           |
|-----------------|-----------------------------------------------------------|
| NBP2-27231PCP   | Mouse IgG2b Isotype Control (MPC-11) [PerCP]              |
| NBP2-59553-50ug | Recombinant Human t-Plasminogen Activator/tPA His Protein |
| 210-TA-005      | TNF-alpha [Unconjugated]                                  |
| DTPA00          | t-Plasminogen Activator/tPA [HRP]                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-74568PCP](http://www.novusbio.com/reviews/submit/NBP2-74568PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

