

# Product Datasheet

## Lumican Antibody (JE11-45)

### NBP2-76847

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

#### Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-76847](http://www.novusbio.com/NBP2-76847)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-76847](http://www.novusbio.com/reviews/destination/NBP2-76847)



**NBP2-76847**

Lumican Antibody (JE11-45)


| Product Information     |                                                                                        |
|-------------------------|----------------------------------------------------------------------------------------|
| Unit Size               | 100 ul                                                                                 |
| Concentration           | 1 mg/ml                                                                                |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality               | Monoclonal                                                                             |
| Clone                   | JE11-45                                                                                |
| Preservative            | 0.05% Sodium Azide                                                                     |
| Isotype                 | IgG                                                                                    |
| Purity                  | Protein A purified                                                                     |
| Buffer                  | TBS (pH7.4), 0.05% BSA, 40% Glycerol                                                   |
| Target Molecular Weight | 38 kDa                                                                                 |

| Product Description |                                                                                                                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit Lumican Antibody (JE11-45) (NBP2-76847) is a recombinant monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. Anti-Lumican Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                | Rabbit                                                                                                                                                                                                                                                                   |
| Gene ID             | 4060                                                                                                                                                                                                                                                                     |
| Gene Symbol         | LUM                                                                                                                                                                                                                                                                      |
| Species             | Human, Mouse, Rat                                                                                                                                                                                                                                                        |
| Immunogen           | Recombinant protein within Human Lumican aa 172-338 / 338. (SwissProt: P51884 Human)                                                                                                                                                                                     |


| Product Application Details |                                                                                                                                                                          |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                               |
| Recommended Dilutions       | Western Blot 1:500-1:2000, Flow Cytometry 1:50-1:100, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin 1:50-1:200 |

**Images**


Western Blot: Lumican Antibody (JE11-45) [NBP2-76847] - Western blot analysis of Lumican on different lysates using anti-Lumican antibody at 1/1,000 dilution. Positive control: Lane 1: Jurkat Lane 2: 293T Lane 3: Human skin




Immunocytochemistry/Immunofluorescence: Lumican Antibody (JE11-45) [NBP2-76847] - ICC staining Lumican in MCF-7 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Immunohistochemistry: Lumican Antibody (JE11-45) [NBP2-76847] - Immunohistochemical analysis of paraffin-embedded human skin tissue using anti-Lumican antibody. Counter stained with hematoxylin.




Flow Cytometry: Lumican Antibody (JE11-45) [NBP2-76847] - Flow cytometric analysis of 293T cells with Lumican antibody at 1/100 dilution (purple) compared with an unlabelled control (cells without incubation with primary antibody; yellow). Alexa Fluor 488-conjugated goat anti-rabbit IgG was used as the secondary antibody.




Immunocytochemistry/Immunofluorescence: Lumican Antibody (JE11-45) [NBP2-76847] - ICC staining Lumican in MG-63 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunocytochemistry/Immunofluorescence: Lumican Antibody (JE11-45) [NBP2-76847] - ICC staining Lumican in PANC-1 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.



Immunohistochemistry: Lumican Antibody (JE11-45) [NBP2-76847] - Immunohistochemical analysis of paraffin-embedded rat kidney tissue using anti-Lumican antibody. Counter stained with hematoxylin.



## Publications

Barnes EA, Ito R, Che X et al. Loss of Prolyl Hydroxylase 1 and 2 in SM22<sup>+</sup>-Expressing Cells Prevents Hypoxia-Induced Pulmonary Hypertension American journal of physiology. Lung cellular and molecular physiology 2023-10-17 [PMID: 37847687]

Alves AL, Carvalho AC, Machado I et al. Cell-Laden Marine Gelatin Methacryloyl Hydrogels Enriched with Ascorbic Acid for Corneal Stroma Regeneration Bioengineering (Basel, Switzerland) 2023-01-04 [PMID: 36671634] (ICC/IF, Human)

Details:  
Dilution used in ICC 1:200



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-76847**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-87726PEP | Lumican Recombinant Protein Antigen                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-76847](http://www.novusbio.com/reviews/submit/NBP2-76847)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

