

# Product Datasheet

## Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2) [HRP] NBP2-80053H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-80053H](http://www.novusbio.com/NBP2-80053H)

Updated 12/3/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-80053H](http://www.novusbio.com/reviews/destination/NBP2-80053H)



**NBP2-80053H**

Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2) [HRP]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 4H8 WH2                                                                                     |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG2a                                                                                       |
| Conjugate           | HRP                                                                                         |
| Purity              | Protein G purified                                                                          |
| Buffer              | PBS                                                                                         |

| Product Description |                                                                                |
|---------------------|--------------------------------------------------------------------------------|
| Host                | Rat                                                                            |
| Gene ID             | 4055                                                                           |
| Gene Symbol         | LTBR                                                                           |
| Species             | Mouse                                                                          |
| Immunogen           | Recombinant mouse Lymphotoxin beta R/TNFRSF3 (cysteine-rich region/aa 31-221). |

| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Flow Cytometry, Functional Assay                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Recommended Dilutions       | Flow Cytometry, Functional Assay                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Application Notes           | The monoclonal antibody to mouse LT-beta-R is an agonist that can be used for the investigation of the regulation of BAFF (BlyS), chemokines and integrins using in vivo and tissue culture models, the development of NK cells and NK T cells, to study the regulation of NF-kappa-B family of transcription factors in regulation of inflammation and homeostasis, particularly RelB NF-kappa-B2 pathway. For use as an agonist the MAb to LT-beta-R is added to cell cultures at 2ug/ml. For in vivo use, mice are injected intraperitoneally with 50ug of agonistic MAb to LT-beta-R in sterile phosphate saline buffer. Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-80053H**

---

|                   |                                                          |
|-------------------|----------------------------------------------------------|
| NBP2-51906-0.05mg | Recombinant Human Lymphotoxin beta R/TNFRSF3 His Protein |
| 210-TA-005        | TNF-alpha [Unconjugated]                                 |
| 7538-LR-100       | Lymphotoxin beta R/TNFRSF3 [Unconjugated]                |
| D6050             | IL-6 [HRP]                                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-80053H](http://www.novusbio.com/reviews/submit/NBP2-80053H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

