

Product Datasheet

Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2) - Azide and BSA Free NBP2-80054

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-80054

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-80054



NBP2-80054

Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2) - Azide and BSA Free

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4H8 WH2
Preservative	No Preservative
Isotype	IgG2a
Purity	Protein G purified
Buffer	PBS containing 10% glycerol
Product Description	
Description	Novus Biologicals Rat Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2) - Azide and BSA Free (NBP2-80053) is a monoclonal antibody validated for use in Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rat
Gene ID	4055
Gene Symbol	LTBR
Species	Mouse
Immunogen	Recombinant mouse Lymphotoxin beta R/TNFRSF3 (cysteine-rich region/aa 31-221).
Product Application Details	
Applications	Flow Cytometry, Functional Assay
Recommended Dilutions	Flow Cytometry, Functional Assay
Application Notes	The monoclonal antibody to mouse LT-beta-R is an agonist that can be used for the investigation of the regulation of BAFF (BlyS), chemokines and integrins using in vivo and tissue culture models, the development of NK cells and NK T cells, to study the regulation of NF-kappa-B family of transcription factors in regulation of inflammation and homeostasis, particularly RelB NF-kappa-B2 pathway. For use as an agonist the MAb to LT-beta-R is added to cell cultures at 2ug/ml. For in vivo use, mice are injected intraperitoneally with 50ug of agonistic MAb to LT-beta-R in sterile phosphate saline buffer.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-80054

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
F0105B	Goat anti-Rat IgG Secondary Antibody [Phycoerythrin]
NBP2-21947-0.1mg	Rat IgG2a Isotype Control (2A3)
NBP2-51906-0.05mg	Recombinant Human Lymphotoxin beta R/TNFRSF3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-80054

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

