

# Product Datasheet

## **Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2)** **[Janelia Fluor® 646]** **NBP2-80054JF646**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-80054JF646](http://www.novusbio.com/NBP2-80054JF646)

Updated 12/3/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-80054JF646](http://www.novusbio.com/reviews/destination/NBP2-80054JF646)



**NBP2-80054JF646**

Lymphotoxin beta R/TNFRSF3 Antibody (4H8 WH2) [Janelia Fluor® 646]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	4H8 WH2
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	Janelia Fluor 646
Purity	>95%, by SDS-PAGE
Buffer	50mM Sodium Borate
Product Description	
Host	Rat
Gene ID	4055
Gene Symbol	LTBR
Species	Mouse
Immunogen	Recombinant mouse Lymphotoxin beta R/TNFRSF3 (cysteine-rich region/aa 31-221).
Notes	Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.
Product Application Details	
Applications	Flow Cytometry, Functional Assay
Recommended Dilutions	Flow Cytometry, Functional Assay
Application Notes	The monoclonal antibody to mouse LT-beta-R is an agonist that can be used for the investigation of the regulation of BAFF (BlyS), chemokines and integrins using in vivo and tissue culture models, the development of NK cells and NK T cells, to study the regulation of NF-kappa-B family of transcription factors in regulation of inflammation and homeostasis, particularly RelB NF-kappa-B2 pathway. For use as an agonist the MAb to LT-beta-R is added to cell cultures at 2ug/ml. For in vivo use, mice are injected intraperitoneally with 50ug of agonistic MAb to LT-beta-R in sterile phosphate saline buffer. Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-80054JF646**

---

NBP2-51906-0.05mg	Recombinant Human Lymphotoxin beta R/TNFRSF3 His Protein
210-TA-005	TNF-alpha [Unconjugated]
7538-LR-100	Lymphotoxin beta R/TNFRSF3 [Unconjugated]
D6050	IL-6 [HRP]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-80054JF646](http://www.novusbio.com/reviews/submit/NBP2-80054JF646)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

