

Product Datasheet

CD27/TNFRSF7 Antibody (LG.3A10) - Azide and BSA Free NBP2-81095-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81095

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81095



NBP2-81095-0.2mg

CD27/TNFRSF7 Antibody (LG.3A10) - Azide and BSA Free

Product Information	
Unit Size	0.2 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	LG.3A10
Preservative	0.02% Proclin 300
Isotype	IgG Kappa
Purity	Protein A purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Armenian Hamster CD27/TNFRSF7 Antibody (LG.3A10) - Azide and BSA Free (NBP2-81095) is a recombinant monoclonal antibody validated for use in Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Armenian Hamster
Gene ID	939
Gene Symbol	CD27
Species	Human, Mouse, Rat
Specificity/Sensitivity	This antibody is specific for mouse CD27/TNFRSF7, which, as in the human system, is largely restricted to cells of the T cell lineage, and is expressed on the great majority of both alpha-beta and gamma-delta T lymphocytes, and thymocytes.
Immunogen	This antibody was produced by injecting Armenian hamsters intraperitoneally with a mixture of two murine CD27/TNFRSF7 tag-expressing Armenian hamster fibroblast clones (AR0524 and AR0530) in Freund's adjuvant

Product Application Details	
Applications	Flow Cytometry, Functional, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Functional
Application Notes	This antibody is specific for murine CD27, and has been used to assess the distribution of CD27 expression using IF, FC and IP analysis (Gravstein, 1995). Upon cross-linking, this antibody has been shown to induce a 4-fold increase in T cell proliferation in response to suboptimal stimulation with concanavalin A (Gravstein, 1995), indicating that this antibody is able to mimic ligand binding.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81095-0.2mg

NB120-5738	Goat anti-Armenian Hamster IgG (H+L) Secondary Antibody
NBP2-38434PEP	CD27/TNFRSF7 Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
DY382-05	CD27/TNFRSF7 [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81095

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

