

Product Datasheet

B7-1/CD80 Antibody (IDEC-114 (Galiximab)) - Chimeric - Azide and BSA Free NBP2-81100-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81100

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81100



NBP2-81100-0.2mg

B7-1/CD80 Antibody (IDEC-114 (Galiximab)) - Chimeric - Azide and BSA Free

| Product Information | |
|---------------------|--|
| Unit Size | 0.2 mg |
| Concentration | 1 mg/ml |
| Storage | Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C. |
| Clonality | Monoclonal |
| Clone | IDEC-114 (Galiximab) |
| Preservative | 0.02% Proclin 300 |
| Isotype | IgG Lambda |
| Purity | Protein A purified |
| Buffer | PBS |


| Product Description | |
|-------------------------|--|
| Description | Novus Biologicals Rabbit B7-1/CD80 Antibody (IDEC-114 (Galiximab)) - Chimeric - Azide and BSA Free (NBP2-81100) is a recombinant monoclonal antibody validated for use in IHC, WB, ELISA, Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 941 |
| Gene Symbol | CD80 |
| Species | Human |
| Specificity/Sensitivity | Galiximab is a primatized mAb which consists of human constant and primate (cynomolgus macaque) variable regions and binds specifically to B7-1/CD80. B7-1/CD80 is a surface glycoprotein and a member of the B7 family of costimulatory molecules. B7-1/CD80 antigen regulates T cell activation (through interacting with CD28 or CD152) and is expressed transiently in antigen-presenting cells, T cells and normal B cells, and expressed constitutively on various subtypes of B-cell lymphomas. |
| Immunogen | Galiximab was prepared by immunizing cynomolgus monkeys with recombinant B7-1/CD80 antigen. The variable regions of the light and heavy chains were then cloned by being incorporated into a cassette vector (N5LG1) containing human constant region genes and subsequently transfected into the Dg44 CHO cell line (Hariharan K, 2013) |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, ELISA, Flow Cytometry, Functional, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Block/Neutralize |
| Recommended Dilutions | Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Functional, Block/Neutralize |
| Application Notes | This chimeric rabbit antibody was made using the variable domain sequences of the original Human IgG1 format, for improved compatibility with existing reagents, assays and techniques. Galiximab acts to decrease in cell proliferation, inhibition of the constitutively active NF-kappaB pathway, increase in apoptosis and ADCC against various B-cell lymphoma cell lines. Galiximab affects CD28 and CD152 interactions with CD80. This antibody can also be used for FC, WB, ELISA, IP and IHC. |



Images

Immunocytochemistry/Immunofluorescence: B7-1/CD80 Antibody (IDEC-114 (Galiximab)) - Chimeric [NBP2-81100] - Immunofluorescence staining of fixed Daudi cells with B7-1/CD80 Antibody (IDEC-114 (Galiximab)) [NBP2-81100]. Immunofluorescence analysis of paraformaldehyde fixed Daudi cells on Shi-fix coverslips, permeabilized with 0.15% Triton and stained with the chimeric rabbit IgG version of IDEC-114 (Galiximab) [NBP2-81100] at 10 ug/ml for 1h followed by Alexa Fluor 488 secondary antibody (1 ug/ml), showing membrane staining. The nuclear stain is DAPI (blue). Panels show from left-right, top-bottom NBP2-81100, DAPI, merged channels and an isotype control. The isotype control was stained with an anti-Fluorescein antibody followed by Alexa Fluor 488 secondary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81100-0.2mg

| | |
|------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-26578 | Recombinant Mouse B7-1/CD80 Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81100

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

