

Product Datasheet

Human INSL3 ELISA Kit (Colorimetric) NBP2-81140

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-81140

Updated 9/23/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-81140



NBP2-81140**Human INSL3 ELISA Kit (Colorimetric)****Product Information**

Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.

Product Description


Gene ID	3640
Gene Symbol	INSL3
Species	Human
Specificity/Sensitivity	This kit recognizes Human INSL3 in samples. No significant cross-reactivity or interference between Human INSL3 and analogues was observed.
Kit Components	Micro ELISA Plate(Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Reference Standard & Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Description, Certificate of Analysis
Standard Curve Range	0.16 - 10 ng/mL
Sensitivity	0.09 ng/mL
Inter Assay	CV% < 5.48%
Intra Assay	CV% < 5.18%
Assay Type	Sandwich-ELISA
Spiking Recovery	92-109%
Suitable Sample Type	Serum, plasma, cell supernatant, and other biological fluids.
Sample Volume	100 ul

Product Application Details

Applications	ELISA
Recommended Dilutions	ELISA

Images

Human INSL3 ELISA Kit (Colorimetric) [NBP2-81140] - Typical standard curve. The provided standard curve is for reference only, as the OD values may vary according to the conditions of the actual experiment. The standard curve should be established for each test.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-81140

NBP1-81223PEP	INSL3 Recombinant Protein Antigen
202-IL-010	IL-2 [Unconjugated]
NBP1-81223	INSL3 Antibody
DY262	Amphiregulin [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-81140

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

