

Product Datasheet

Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 Antibody [DyLight 350] NBP2-82107UV

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-82107UV

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-82107UV



NBP2-82107UV

Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 Antibody [DyLight 350]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 350
Purity	Peptide affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	23135
Gene Symbol	KDM6B
Species	Human, Mouse, Rat
Reactivity Notes	0
Specificity/Sensitivity	This antibody is specific for Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 and will not recognize other JMJD proteins. Human Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 has 2 isoforms, including isoform 1 (1682aa, 181kD) and isoform 2 (1643aa, 177kD). Mouse Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 has one isoform (1641aa, 176kD) and rat Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 also has one isoform (1356aa, 146kD). NBP2-82107 can detect human, mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3. The immunogen is located within amino acids 1370-1420 of Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3. Amino Acid Sequence: KVPGSRTPGHQENNNF
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunocytochemistry
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunocytochemistry
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-82107UV

NBP2-24891UV	Rabbit IgG Isotype Control [DyLight 350]
NBP1-06640PEP	Lysine (K)-specific Demethylase 6B/KDM6B/JMJD3 Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
NB200-106	p19ARF/CDKN2A Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-82107UV

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

