

# Product Datasheet

## LRRTM3 Antibody - BSA Free

### NBP2-82277-0.1ml

Unit Size: 0.1ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-82277](http://www.novusbio.com/NBP2-82277)

Updated 2/23/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-82277](http://www.novusbio.com/reviews/destination/NBP2-82277)



**NBP2-82277-0.1ml**

LRRTM3 Antibody - BSA Free

**Product Information**

|                                |  |
|--------------------------------|--|
| <b>Unit Size</b>               | 0.1ml  |
| <b>Concentration</b>           | 0.5 mg/ml  |
| <b>Storage</b>                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>               | Polyclonal   |
| <b>Preservative</b>            | 0.09% Sodium Azide   |
| <b>Purity</b>                  | Affinity purified  |
| <b>Buffer</b>                  | PBS, 2% Sucrose  |
| <b>Target Molecular Weight</b> | 56 kDa   |

**Product Description**


|                    |  |
|--------------------|--|
| <b>Host</b>        | Rabbit   |
| <b>Gene ID</b>     | 347731   |
| <b>Gene Symbol</b> | LRRTM3   |
| <b>Species</b>     | Human  |
| <b>Immunogen</b>   | The immunogen is a synthetic peptide directed towards the middle region of Human LRRTM3. Peptide sequence:<br>QLQQRSLMRRHRKKKRQSLKQMTPTSTQEFYVDYKPTNTETSEMLLNGTGP<br>The peptide sequence for this immunogen was taken from within the described region. |

**Product Application Details**

|                              |                        |
|------------------------------|------------------------|
| <b>Applications</b>          | Western Blot           |
| <b>Recommended Dilutions</b> | Western Blot 1.0 ug/ml |

**Images**

Western Blot: LRRTM3 Antibody [NBP2-82277] - Host: Rabbit. Target Name: LRRTM3. Sample Type: HT1080 Whole Cell lysates. Antibody Dilution: 1.0ug/ml





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-82277-0.1ml**

---

|               |   |
|---------------|---|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP3-21335PEP | LRRTM3 Recombinant Protein Antigen                  |
| 1310-SE-010   | u-Plasminogen Activator/Urokinase [Unconjugated]    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-82277](http://www.novusbio.com/reviews/submit/NBP2-82277)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

