

Product Datasheet

Methyltransferase like 3 Antibody - BSA Free NBP2-82284-0.1ml

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-82284

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-82284



NBP2-82284-0.1ml

Methyltransferase like 3 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Protein A purified
Buffer	PBS, 2% Sucrose
Target Molecular Weight	64 kDa


Product Description	
Description	Novus Biologicals Rabbit Methyltransferase like 3 Antibody - BSA Free (NBP2-82284) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	56339
Gene Symbol	METTL3
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human Methyltransferase like 3. Peptide sequence: MSDTWSSIQAHHKQLDSLRLQRRRKQDSGHLDLRNPEAALSPTFRSDS The peptide sequence for this immunogen was taken from within the described region.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin




Images

Western Blot: Methyltransferase like 3 Antibody [NBP2-82284] - WB
Suggested Anti-METTL3 Antibody Titration: 5.0ug/ml. ELISA Titer:
1:312500. Positive Control: Human Thymus



Immunohistochemistry: Methyltransferase like 3 Antibody [NBP2-82284]
- Human Intestine



Western Blot: Methyltransferase like 3 Antibody [NBP2-82284] - Host:
Rabbit. Target: METTL3. Positive control (+): HepG2 Cell Lysate (HG).
Negative control (-): Human Liver Tumor (T-LI). Antibody concentration:
1ug/ml



Host: Rabbit
Target name: METTL3
Positive control: ~25ug HepG2 Cell Lysate (HG)
Negative control: ~25ug Human Liver (LI)
Antibody concentration: 1 ug/ml



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-82284-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
H00056339-P01-10ug	Recombinant Human Methyltransferase like 3 GST (N-Term) Protein
DBD00	BDNF [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-82284

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

