

Product Datasheet

FAM108A1 Antibody - BSA Free NBP2-84877-0.1ml

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-84877

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-84877



NBP2-84877-0.1ml

FAM108A1 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose

Product Description


Description	Novus Biologicals Rabbit FAM108A1 Antibody - BSA Free (NBP2-84877) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	81926
Gene Symbol	ABHD17A
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human FAM108A1. Peptide sequence: VPTVDLASRYECAAVVLHSPLTSGMRVAFPDTKKTYCFDAFPNIEKVSKI The peptide sequence for this immunogen was taken from within the described region.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml

Images

Western Blot: FAM108A1 Antibody [NBP2-84877] - Host: Rabbit. Target Name: FAM108A1. Sample Type: Human Adult Placenta. Antibody Dilution: 1.0ug/ml


FAM108A1

Western Blot: FAM108A1 Antibody [NBP2-84877] - Host: Rabbit. Target Name: FAM108A1. Sample Type: Fetal Brain lysates. Antibody Dilution: 1.0ug/ml




Western Blot: FAM108A1 Antibody [NBP2-84877] - Host: Rabbit. Target Name: FAM108A1. Sample Type: Human Fetal Lung. Antibody Dilution: 1.0ug/ml

FAM108A1




Western Blot: FAM108A1 Antibody [NBP2-84877] - Host: Rabbit. Target Name: FAM108A1. Sample Type: Human Fetal Heart. Antibody Dilution: 1.0ug/ml

FAM108A1




Western Blot: FAM108A1 Antibody [NBP2-84877] - Host: Rabbit. Target Name: FAM108A1. Sample Type: Human 293T. Antibody Dilution: 1.0ug/ml
ABHD17A is supported by BioGPS gene expression data to be expressed in HEK293T

FAM108A1



Western Blot: FAM108A1 Antibody [NBP2-84877] - Host: Rabbit. Target Name: FAM108A1. Sample Type: Human Jurkat. Antibody Dilution: 1.0ug/ml ABHD17A is supported by BioGPS gene expression data to be expressed in Jurkat

FAM108A1



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-84877-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
H00081926-P01-10ug	Recombinant Human FAM108A1 GST (N-Term) Protein
5255-EN-010	ENPP-2/Autotaxin [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-84877

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

