

# Product Datasheet

## PITX3 Antibody - BSA Free

### NBP2-85479-0.1ml

Unit Size: 0.1ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-85479](http://www.novusbio.com/NBP2-85479)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-85479](http://www.novusbio.com/reviews/destination/NBP2-85479)



**NBP2-85479-0.1ml**


PITX3 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1ml                                                                                                                                                                                                                                           |
| Concentration               | 0.5 mg/ml                                                                                                                                                                                                                                       |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                                                                                          |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                      |
| Preservative                | 0.09% Sodium Azide                                                                                                                                                                                                                              |
| Purity                      | Affinity purified                                                                                                                                                                                                                               |
| Buffer                      | PBS, 2% Sucrose                                                                                                                                                                                                                                 |
| Product Description         |                                                                                                                                                                                                                                                 |
| Description                 | Novus Biologicals Rabbit PITX3 Antibody - BSA Free (NBP2-85479) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                           |
| Host                        | Rabbit                                                                                                                                                                                                                                          |
| Gene ID                     | 5309                                                                                                                                                                                                                                            |
| Gene Symbol                 | PITX3                                                                                                                                                                                                                                           |
| Species                     | Human                                                                                                                                                                                                                                           |
| Immunogen                   | The immunogen is a synthetic peptide directed towards the middle region of human PITX3. Peptide sequence: YEEVYPGYSYGNWPPKALAPPLAAKTFPFAFNSVNVGPLASQPVFSPSS The peptide sequence for this immunogen was taken from within the described region. |
| Product Application Details |                                                                                                                                                                                                                                                 |
| Applications                | Western Blot                                                                                                                                                                                                                                    |
| Recommended Dilutions       | Western Blot 1.0 ug/ml                                                                                                                                                                                                                          |



## Images


Western Blot: PITX3 Antibody [NBP2-85479] - Host: Rabbit. Target Name: PITX3. Sample Type: Human Fetal Liver. Antibody Dilution: 1.0ug/ml



Western Blot: PITX3 Antibody [NBP2-85479] - WB Suggested Anti-PITX3 Antibody Titration: 0.2-1 ug/ml. ELISA Titer: 1:312500. Positive Control: Human Liver



Western Blot: PITX3 Antibody [NBP2-85479] - Host: Rabbit. Target Name: PITX3. Sample Tissue: Human RPMI 8226 Whole Cell. Antibody Dilution: 1ug/ml





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-85479-0.1ml**

---

|            |                                                     |
|------------|-----------------------------------------------------|
| HAF008     | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160     | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| 210-TA-005 | TNF-alpha [Unconjugated]                            |
| NBL1-14449 | PITX3 Overexpression Lysate                         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-85479](http://www.novusbio.com/reviews/submit/NBP2-85479)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

