

# Product Datasheet

## TEAD3 Antibody - BSA Free

### NBP2-85894-0.1ml

Unit Size: 0.1ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-85894](http://www.novusbio.com/NBP2-85894)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-85894](http://www.novusbio.com/reviews/destination/NBP2-85894)



**NBP2-85894-0.1ml**

TEAD3 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                                                  |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1ml                                                                                                                                                                                                                                            |
| Concentration               | 0.5 mg/ml                                                                                                                                                                                                                                        |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                                                                                           |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                       |
| Preservative                | 0.09% Sodium Azide                                                                                                                                                                                                                               |
| Purity                      | Affinity purified                                                                                                                                                                                                                                |
| Buffer                      | PBS, 2% Sucrose                                                                                                                                                                                                                                  |
| Product Description         |                                                                                                                                                                                                                                                  |
| Description                 | Novus Biologicals Rabbit TEAD3 Antibody - BSA Free (NBP2-85894) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                            |
| Host                        | Rabbit                                                                                                                                                                                                                                           |
| Gene ID                     | 7005                                                                                                                                                                                                                                             |
| Gene Symbol                 | TEAD3                                                                                                                                                                                                                                            |
| Species                     | Human                                                                                                                                                                                                                                            |
| Immunogen                   | The immunogen is a synthetic peptide directed towards the middle region of human TEAD3. Peptide sequence: KFSPPSPLPQAVFSTSSRFWSSPPLLGGQPGPSQDIKPFAQPAYPIQPPL The peptide sequence for this immunogen was taken from within the described region. |
| Product Application Details |                                                                                                                                                                                                                                                  |
| Applications                | Western Blot                                                                                                                                                                                                                                     |
| Recommended Dilutions       | Western Blot 1.0 ug/ml                                                                                                                                                                                                                           |




## Images

Western Blot: TEAD3 Antibody [NBP2-85894] - WB Suggested Anti-TEAD3 Antibody Titration: 0.2-1 ug/ml. ELISA Titer: 1:62500. Positive Control: Transfected 293T



Western Blot: TEAD3 Antibody [NBP2-85894] - WB Suggested Anti-TEAD3 antibody Titration: 1 ug/mL. Sample Type: Human heart





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-85894-0.1ml**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| H00007005-P01-10ug | Recombinant Human TEAD3 GST (N-Term) Protein        |
| NBP1-83949PEP      | TEAD3 Recombinant Protein Antigen                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-85894](http://www.novusbio.com/reviews/submit/NBP2-85894)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

