

Product Datasheet

ADA2a Antibody - BSA Free

NBP2-86955

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-86955

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-86955



NBP2-86955


ADA2a Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Product Description	
Description	Novus Biologicals Rabbit ADA2a Antibody - BSA Free (NBP2-86955) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6871
Gene Symbol	TADA2A
Species	Human
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of human ADA2a. Peptide sequence: MDRLGPFSNDPSDKPPCRGCGSSYLMEPYIKCAECGPPPPFFLCLQCFTTRGF The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml




Images

Western Blot: ADA2a Antibody [NBP2-86955] - WB Suggested Anti-TADA2L Antibody Titration: 0.2-1 ug/ml. ELISA Titer: 1:312500. Positive Control: Jurkat cell lysate. TADA2A is strongly supported by BioGPS gene expression data to be expressed in Jurkat



Western Blot: ADA2a Antibody [NBP2-86955] - Host: Rabbit. Target Name: TADA2A. Sample Type: 721_B. Antibody Dilution: 1.0ug/ml. TADA2A is strongly supported by BioGPS gene expression data to be expressed in Human 721_B cells

TADA2A





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-86955

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-57507PEP	ADA2a Recombinant Protein Antigen
2695-SE-010	Airway Trypsin-like Protease/HAT/TMPRSS11D

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-86955

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

