

# Product Datasheet

## HES6 Antibody - BSA Free NBP2-87563

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-87563](http://www.novusbio.com/NBP2-87563)

Updated 1/6/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-87563](http://www.novusbio.com/reviews/destination/NBP2-87563)



**NBP2-87563**


HES6 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                                                        |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                                                                                 |
| Concentration               | 1 mg/ml                                                                                                                                                                                                                                                |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                                                                                                 |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                             |
| Preservative                | 0.09% Sodium Azide                                                                                                                                                                                                                                     |
| Purity                      | Protein A purified                                                                                                                                                                                                                                     |
| Buffer                      | PBS, 2% Sucrose                                                                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                                                                        |
| Description                 | Novus Biologicals Rabbit HES6 Antibody - BSA Free (NBP2-87563) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                           |
| Host                        | Rabbit                                                                                                                                                                                                                                                 |
| Gene ID                     | 55502                                                                                                                                                                                                                                                  |
| Gene Symbol                 | HES6                                                                                                                                                                                                                                                   |
| Species                     | Human                                                                                                                                                                                                                                                  |
| Immunogen                   | The immunogen is a synthetic peptide directed towards the C-terminal region of human HES6. Peptide sequence: PPGPGDDLCSDL EEAP EAELSQA PAEGPDLVPAALGSLTTAQIARSVWRP The peptide sequence for this immunogen was taken from within the described region. |
| Product Application Details |                                                                                                                                                                                                                                                        |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot 1.0 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                            |




## Images

Western Blot: HES6 Antibody [NBP2-87563] - WB Suggested Anti-HES6  
Antibody Titration: 1.25ug/ml. Positive Control: Transfected 293T



Immunohistochemistry: HES6 Antibody [NBP2-87563] - Human  
Pancreas





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-87563](http://www.novusbio.com/reviews/submit/NBP2-87563)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

