

Product Datasheet

MSH6 Antibody (035) - Azide and BSA Free NBP2-89265-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-89265

Updated 2/26/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-89265



NBP2-89265-100ul

MSH6 Antibody (035) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	035
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	0.2 um filtered solution in PBS

Product Description	
Description	This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Gene ID	2956
Gene Symbol	MSH6
Species	Human, Mouse, Rat, Ferret, Monkey
Reactivity Notes	Cross reactivity in IHC: Cynomolgus, Mouse, Rat, Ferret
Immunogen	This antibody was obtained from a rabbit immunized with a synthetic peptide corresponding to the center region of the Human MSH6.


Product Application Details	
Applications	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:100-1:500, Immunocytochemistry/ Immunofluorescence 1:20-1:100, Immunohistochemistry-Paraffin 1:100-1:500

Images


Immunocytochemistry/Immunofluorescence: MSH6 Antibody (035) [NBP2-89265] - Staining of MSH6 in HeLa cells. Cells were fixed with 4% PFA, permeabilized with 0.1% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-human MSH6 monoclonal antibody (dilution ratio 1:60) at 4c overnight. Then cells were stained with the Alexa Fluor(R)488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Nucleus.




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in Ferret kidney with rabbit monoclonal antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in human rectal cancer with rabbit monoclonal antibody (1:200).




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in human colon carcinoma with rabbit monoclonal antibody (1:200).




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in cynomolgus spleen with rabbit monoclonal antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in cynomolgus kidney with rabbit monoclonal antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in mouse spleen with rabbit monoclonal antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in mouse kidney with rabbit monoclonal antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in rat spleen with rabbit monoclonal antibody at 1:200 dilution.



Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in rat kidney with rabbit monoclonal antibody at 1:200 dilution



Immunohistochemistry-Paraffin: MSH6 Antibody (035) [NBP2-89265] - Staining MSH6 in Ferret spleen with rabbit monoclonal antibody at 1:200 dilution.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-89265-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00002956-Q01-10ug	Recombinant Human MSH6 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-89265

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

