

# Product Datasheet

## Leukotriene A4 Hydrolase/LTA4H Antibody (455) [FITC] NBP2-89375F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-89375F](http://www.novusbio.com/NBP2-89375F)

Updated 11/12/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-89375F](http://www.novusbio.com/reviews/destination/NBP2-89375F)



**NBP2-89375F**

Leukotriene A4 Hydrolase/LTA4H Antibody (455) [FITC]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | 455   |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG   |
| Conjugate           | FITC  |
| Purity              | Protein A purified  |
| Buffer              | PBS   |

| Product Description     |  |
|-------------------------|--|
| Host                    | Rabbit   |
| Gene ID                 | 4048   |
| Gene Symbol             | LTA4H  |
| Species                 | Human  |
| Specificity/Sensitivity | Has cross-reactivity in ELISA with:<br><br>Mouse LTA4H   |
| Immunogen               | This antibody was obtained from a rabbit immunized with purified, recombinant Human Leukotriene A4 Hydrolase/LTA4H (Accession#: NP_000886.1; Met 1-Asp 611). |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | ELISA  |
| Recommended Dilutions       | ELISA  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-89375F**

---

|               |  |
|---------------|--|
| NBP2-24892    | Rabbit IgG Isotype Control [FITC]                          |
| NBP1-90375PEP | Leukotriene A4 Hydrolase/LTA4H Recombinant Protein Antigen |
| 210-TA-005    | TNF-alpha [Unconjugated]                                   |
| 4008-ZN-010   | Leukotriene A4 Hydrolase/LTA4H [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-89375F](http://www.novusbio.com/reviews/submit/NBP2-89375F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

