

Product Datasheet

Serpin A8/Angiotensinogen Antibody (1) [mFluor Violet 610 SE] NBP2-89736MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-89736MFV610

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-89736MFV610



NBP2-89736MFV610

Serpine A8/Angiotensinogen Antibody (1) [mFluor Violet 610 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	1
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	mFluor Violet 610 SE
Purity	Protein A purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	183
Gene Symbol	AGT
Species	Human
Specificity/Sensitivity	No cross-reactivity in ELISA with: Human SerpinA1 / A1AT Human SerpinA3 / SERPINA3 / AACT Human SerpinA5 / Protein C inhibitor Human SerpinD1 / HC2 Human SerpinF2 / SERPINF2 Human SerpinG1 / C1IN Human SerpinI1 / Protease inhibitor 12 Human SerpinC1 /SERPINC1/ Antithrombin-III
Immunogen	This antibody was obtained from a rabbit immunized with purified, recombinant Human Serpin A8/Angiotensinogen (Accession#: NP_000020.1; Met 1-Ala 485).
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Images

Serpin A8/Angiotensinogen Antibody (1) [mFluor Violet 610 SE] [NBP2-89736MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-89736MFV610

NBP1-30027PEP	Serpin A8/Angiotensinogen Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
DY3156-05	Serpin A8/Angiotensinogen [Biotin]
NBP2-22203	ERK1 Antibody (1E5)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-89736MFV610

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

