

Product Datasheet

DNTTIP2 Antibody - Azide and BSA Free NBP2-92296-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-92296

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-92296



NBP2-92296-0.1ml


DNTTIP2 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	84 kDa
Product Description	
Description	Novus Biologicals Rabbit DNTTIP2 Antibody - Azide and BSA Free (NBP2-92296) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	30836
Gene Symbol	DNTTIP2
Species	Human, Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human DNTTIP2 (NP_055412.2). M V V T R S A R A K A S I Q A A S A E S S G Q K S F A A N G I Q A H P E S S T G S D A R T T A E S Q T T G K Q S L I P R T P K A R K R K S R T T G S L P K G T E P S T D G E T S E A E S N Y S V S E H H D T I L R V T R R R Q I L I A C S P V S S V R K K P K V T P T K E S Y T E E I V S E A E S H V S G I S R I V L P T E K T T G A R R S K A K S L T D P S Q E S H T E A I S D A E T S S S D I S F S G I A T R R T R S M Q R K L K A Q T E K K D S K I V P G N E K Q I V G T P V N S E D S D T R Q T S H L Q A R S L S E I N K P N F Y N N D F D D D F S H R S S E N I L T V H E Q A N V E S L K E T K Q N C K D L D
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500-1:2000, Immunocytochemistry/ Immunofluorescence 1:50-1:200




Images

Western Blot: DNTTIP2 Antibody [NBP2-92296] - Analysis of extracts of mouse pancreas, using DNTTIP2 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 10s.



Immunocytochemistry/Immunofluorescence: DNTTIP2 Antibody [NBP2-92296] - Analysis of U-2 OS cells using DNTTIP2 at dilution of 1:100. Blue: DAPI for nuclear staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-92296-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-38710PEP	DNTTIP2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-92296

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

