

Product Datasheet

CEP164 Antibody - BSA Free

NBP2-92844-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-92844

Updated 4/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-92844



NBP2-92844-0.1ml

CEP164 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol

Product Description


Host	Rabbit
Gene ID	22897
Gene Symbol	CEP164
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 730-829 of human CEP164 (NP_055771.4). EASEKSEQAALNAAKEKALQQLREQLEGERKEAVATLEKEHSAELERLCSSLE AKHREVVSSLQKKIQEAQQKEEAQLQKCLGQVEHRVHQKSYHVAGYE

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:2000, Immunohistochemistry 1:50-1:200, Immunohistochemistry-Paraffin

Images


Western Blot: CEP164 Antibody [NBP2-92844] - Analysis of extracts of A-375 cells, using CEP164 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time




Immunohistochemistry-Paraffin: CEP164 Antibody [NBP2-92844] -
Paraffin-embedded rat liver using CEP164 .




Immunohistochemistry-Paraffin: CEP164 Antibody [NBP2-92844] -
Paraffin-embedded human liver damage using CEP164 .



Immunohistochemistry-Paraffin: CEP164 Antibody [NBP2-92844] -
Paraffin-embedded mouse brain using CEP164 .



Immunohistochemistry-Paraffin: CEP164 Antibody [NBP2-92844] -
Paraffin-embedded rat heart using CEP164 .





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-92844

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

