

Product Datasheet

POLD3 Antibody - BSA Free

NBP2-93032-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-93032

Updated 11/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-93032



NBP2-93032-0.1ml

POLD3 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol

Product Description

Host	Rabbit
Gene ID	10714
Gene Symbol	POLD3
Species	Human, Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 197-466 of human POLD3 (NP_006582.1). ASKAAAKTQETNKETKTEAKEVTNASAAGNKAPGKGNNMMSNFFGKAAMNKF VNLDSEQAVKEEKIVEQPTVSVTEPKLATPAGLKKSSKKAEPVKVLQKEKKRGK RVALSDDETKETENMRKKRRRIKLPESDSSEDEVFPDSPGAYEAESPSPPPPP SPPLEPVPKTEPEPPSVKSSSGENKRKRKRVLKSPTYLDGEGCIVTEKVYES CTDSEELNMKTSSVHRPPAMTVKKEPREERKGPKKGTAALGKANRQVSITGF FQRK


Product Application Details

Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500 - 1:2000, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/ Immunofluorescence 1:50 - 1:100, Immunohistochemistry-Paraffin 1:100 - 1:500




Images


Western Blot: POLD3 Antibody [NBP2-93032] - Analysis of extracts of various cell lines, using POLD3 . Exposure time: 60s.



Immunocytochemistry/Immunofluorescence: POLD3 Antibody [NBP2-93032] - Analysis of U2OS cells using POLD3 .



Immunocytochemistry/Immunofluorescence: POLD3 Antibody [NBP2-93032] - Analysis of GFP-RNF168 transgenic U2OS cells using POLD3 can be used to mark cells damaged by UV-A laser for they always gather around DNA damage region.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-93032-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00010714-Q01-10ug	Recombinant Human POLD3 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-93032

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

