

# Product Datasheet

## TPK1 Antibody - BSA Free NBP2-93177-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-93177](http://www.novusbio.com/NBP2-93177)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-93177](http://www.novusbio.com/reviews/destination/NBP2-93177)



**NBP2-93177-0.1ml**

TPK1 Antibody - BSA Free

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Affinity purified                                                                           |
| <b>Buffer</b>        | PBS (pH 7.3), 50% glycerol                                                                  |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Novus Biologicals Rabbit TPK1 Antibody - BSA Free (NBP2-93177) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                                   |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Gene ID</b>     | 27010                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Gene Symbol</b> | TPK1                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Species</b>     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen</b>   | Recombinant fusion protein containing a sequence corresponding to amino acids 1-243 of human TPK1 (NP_071890.2).<br>MEHAFTPLELLSTGNLKYCLVILNQPLDNYFRHLWNKALLRACADGGANRLY<br>DITEGERESFLPEFINGDFDSIRPEVREYYATKGCELISTPDQDHTDFTKCLKML<br>QKKIEEKDLKVDVIVTLGGLAGRFDQIMASVNTLFQATHITPFPIIIQEESLIYLLQ<br>PGKHRLHVDTGMEGDWCGLIPVGQPCMQVTTTGLKWNLNTNDVLAFGTLVSTS<br>NTYDGSGVVTVETDHPLLWTMAIKS |

**Product Application Details**

|                              |                             |
|------------------------------|-----------------------------|
| <b>Applications</b>          | Western Blot                |
| <b>Recommended Dilutions</b> | Western Blot 1:500 - 1:1000 |

**Images**

Western Blot: TPK1 Antibody [NBP2-93177] - Western blot analysis of extracts of various cell lines, using TPK1 antibody (NBP2-93177) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-93177-0.1ml**

|                  |                                                     |
|------------------|-----------------------------------------------------|
| HAF008           | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160           | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891       | Rabbit IgG Isotype Control                          |
| NBP1-99102-100ug | Recombinant Human TPK1 His Protein                  |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-93177](http://www.novusbio.com/reviews/submit/NBP2-93177)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

