

# Product Datasheet

## **Pancreatic Lipase Related Protein 2 Antibody - Azide and BSA Free** **NBP2-93552-0.1ml**

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-93552](http://www.novusbio.com/NBP2-93552)

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-93552](http://www.novusbio.com/reviews/destination/NBP2-93552)



**NBP2-93552-0.1ml**


Pancreatic Lipase Related Protein 2 Antibody - Azide and BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                     |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                        |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                      |
| Preservative                | 0.01% Thimerosal                                                                                                                                                                                                                                                                                                |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                             |
| Purity                      | Affinity purified                                                                                                                                                                                                                                                                                               |
| Buffer                      | PBS (pH 7.3), 50% glycerol                                                                                                                                                                                                                                                                                      |
| Target Molecular Weight     | 52 kDa                                                                                                                                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                                                                                                                                 |
| Description                 | Novus Biologicals Rabbit Pancreatic Lipase Related Protein 2 Antibody - Azide and BSA Free (NBP2-93552) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                   |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                          |
| Gene ID                     | 5408                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | PNLIPRP2                                                                                                                                                                                                                                                                                                        |
| Species                     | Mouse, Rat                                                                                                                                                                                                                                                                                                      |
| Immunogen                   | Recombinant fusion protein containing a sequence corresponding to amino acids 301-470 of human PNLIPRP2 (NP_005387.2).<br>LGYPCASYDEFQESKCFPCPAEGCPKMGHYADQFKGKTSAVEQTFFLNTGES<br>GNFTSWRYKVSVTLSGKEKVNGYIRIALYGSNENSKQYEIFKGSLKPDASHTCA<br>IDVDFNVGKIQKVFLWNKRGINLSEPKLGASQITVQSGEDGTEYNFCSSDTVE<br>ENVLQSLYPC |
| Product Application Details |                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                                                                                                                                                                                                                     |
| Recommended Dilutions       | Western Blot 1:500 - 1:1000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:50 - 1:200                                                                                                                                                                                                         |




## Images


Western Blot: Pancreatic Lipase Related Protein 2 Antibody [NBP2-93552] - Analysis of extracts of rat pancreas, using Pancreatic Lipase Related Protein 2 .Exposure time: 90s.



Immunohistochemistry: Pancreatic Lipase Related Protein 2 Antibody [NBP2-93552] - Analysis of rat pancreas using Pancreatic Lipase Related Protein 2 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunohistochemistry: Pancreatic Lipase Related Protein 2 Antibody [NBP2-93552] - Analysis of mouse pancreas using Pancreatic Lipase Related Protein 2 at dilution of 1:100. Blue: DAPI for nuclear staining.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-93552-0.1ml**

---

|                    |                                                                            |
|--------------------|----------------------------------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]                              |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]                        |
| NBP2-24891         | Rabbit IgG Isotype Control                                                 |
| H00005408-P01-10ug | Recombinant Human Pancreatic Lipase Related Protein 2 GST (N-Term) Protein |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-93552](http://www.novusbio.com/reviews/submit/NBP2-93552)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

