

# Product Datasheet

## UBQLNL Antibody - Azide and BSA Free NBP2-93837-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-93837](http://www.novusbio.com/NBP2-93837)

Updated 4/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-93837](http://www.novusbio.com/reviews/destination/NBP2-93837)



**NBP2-93837-0.1ml**


UBQLNL Antibody - Azide and BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                         |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                            |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Preservative                | 0.01% Thimerosal                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Purity                      | Affinity purified                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Buffer                      | PBS (pH 7.3), 50% glycerol                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene ID                     | 143630                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene Symbol                 | UBQLNL                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Immunogen                   | Recombinant fusion protein containing a sequence corresponding to amino acids 176-475 of human UBQLNL (NP_659490.4).<br>LENPSIQRLLSNMEFMWQFISEHLDTQQLMQQNPEVSRLLLDNSEILLQTLELA<br>RNLAMIQEIMQIQPSQNLEYPLNPQPYLGLETMPGGNNALGQNYADINDQML<br>NSMQDPFGGNPFTALLAGQVLEQVQSSPPPPPSQEQDQLTQHPATRVIN<br>SSGGFSSNTSANDTLNKVNHTSKANTAMISTKGQSHICATRQPAWIPALPSIEL<br>TQQLQEEYKDATVSLSSSRQTLKGDQLSDEQSSSQITGGMMQLLMNNPYLA<br>AQIMLFTSMPQLSEQWRQQLPTFLQQTQISDLLSA |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                                                                                                                                                                                                                                                                                                                                                                                               |
| Recommended Dilutions       | Western Blot 1:500-1:2000, Immunocytochemistry/ Immunofluorescence 1:50-1:200                                                                                                                                                                                                                                                                                                                                                                       |




## Images


Western Blot: UBQLNL Antibody [NBP2-93837] - Analysis of extracts of various cell lines, using UBQLNL at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit . Exposure time: 30s.




Immunocytochemistry/Immunofluorescence: UBQLNL Antibody [NBP2-93837] - Analysis of NIH/3T3 cells using UBQLNL at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunocytochemistry/Immunofluorescence: UBQLNL Antibody [NBP2-93837] - Analysis of C6 cells using UBQLNL at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunocytochemistry/Immunofluorescence: UBQLNL Antibody [NBP2-93837] - Analysis of HeLa cells using UBQLNL at dilution of 1:100. Blue: DAPI for nuclear staining.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-93837](http://www.novusbio.com/reviews/submit/NBP2-93837)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

